



OpTest Equipment Inc.

Hawkesbury, ON K6A 3S3
Canada
Phone: 613-632-5169
Fax: 613-632-3744
Email: sales@optest.com

AutoFeed (optional) for HiRes FQA

- Can process 1 to 99 samples automatically for up to 10 hours of continuous running
- User selectable sequence of individual or batch testing
- Automatic & flexible data transfer during or after the testing sequence
- Integrated with HiRes FQA user interface

Note: This option is only available with the HiRes FQA



The AutoFeed allows users to perform fiber analysis on up to 99 samples automatically.

Test samples are prepared in individual 30ml vials. The identification of samples and results are easily entered, stored and viewed with the help of flexible and user-friendly menus similar to the existing HiRes FQA menus. Batch testing is also available. This means that the user is able to average results of selected combinations of vials. It is also possible to have several batches within the 99 samples.

SPECIFICATIONS:

DIMENSIONS [MM]

AutoFeed: 700L x 635W x 178H
(27") (25") (7")

WEIGHT

AutoFeed: 35 kg (74 lbs)

CONNECTIONS

Power: 240/120 V at 50/60 Hz, 600W
stable within 2%

Water: Provided by the HiRes FQA

Instrument Air: Pressure = 680 kPa (100 PSI)
[moisture and oil free]