



## DIGITAL TISSUE ABSORPTION TESTER

Code. X13.100.100

For determination of the water absorption of tissue samples

Applicable standards : ISO 12625-8, DIN EN ISO 12625-



Sold & serviced by:

**OpTest Equipment Inc.**

[www.optest.com](http://www.optest.com) - [sales@optest.com](mailto:sales@optest.com) - +1-613-632-5169



**PAPER TESTING ASSOCIATION**

# LABORATORY Testing Equipment



## SPECIFICATION :

- Chassis made from stainless steel
- Pneumatic lifting of the pot with distilled water
- Automatic detection of the sample weight and water absorption weight
- Detection of the water absorption time
- Basket with 3 g +/- 1 mg acc. standard
- Touch screen display for starting the test and display the values, statistic, curves
- Fully automatic test sequence
- Integrated weighing cell

## CONFIGURATION :

- Instrument with 3 l pot and touch screen display

## DELIVERY CONTENT :

- Absorption tester with built-on touch screen monitor
- Pot for distilled water with 3 l content
- Basket with 3 g weight acc. standard

## CONNECTIONS :

Electrical : 110- 250 V / 50 or 60 Hz

### WEIGHT AND DIMENSIONS (ONLY MACHINE) :

400 x 500 x 700 mm (WxDxH) . 30 Kg

### WEIGHT AND DIMENSIONS WITH PACKING BOX :

500 x 600 x 900 mm (WxDxH) . 50 Kg

## PTA SPAIN

C/ Gurutzegi Nº 16 Pol. Belartza  
20018 SAN SEBASTIÁN  
(Guipúzcoa) - Spain

## PTA FRANCE

14, rue Evariste Galois / ZAC les Ruieres  
38320 EYBENS - France

## PTA ASIA

50 East Coast Road #02-73 Roxy Square  
428769 SINGAPORE - Singapore

[info@ptagroup.com](mailto:info@ptagroup.com)

[www.ptagroup.com](http://www.ptagroup.com)